

Agrifood Technology offers independent testing to the food and pharmaceutical industry, where we can assist manufacturers to verify their raw materials, processes and finished products. Test results can be used to ensure Food Standards Australia New Zealand (FSANZ) and regulations are met.

We report to the limit of reporting (LOR) based on technology and equipment within our laboratory, another option for clients is to report to the limit of detection (LOD). Please contact us for further details on which allergens and specifications are best suited to your needs.

Agrifood Technology recommends customers to perform cleaning validation as part of their manufacturing process. This step is crucial to determine any presence of allergens. We can provide allergen swabs and rinse solutions to validate your processes.

We offer a range of allergens using ELISA kits. The most common allergens are Gluten, Dairy, Peanuts, Crustacean, Egg and Soy. Below is a full list of allergens we currently offer.

Allergen Test Range:

Gluten	Gluten (Fermented)	Total Gluten	Egg
Total Milk Protein	Beta Lacto Globulin	Casein	Soy
Peanut	Lupin	Pine Nut	Oat Protein
Crustacean	Fish	Mollusc	Sesame
TREE NUTS:			
Almond	Brazil Nut	Cashew	Hazelnut
Macadamia	Pecan	Pistachio	Walnut

Agrifood Technology recommends regular allergen testing to minimise risk and to provide our customers with confidence and ease of mind.

CONTACT US:

VIC Laboratory: lab.vic@agrifood.com.au / 1800 801 312 WA Laboratory: lab.wa@agrifood.com.au / (08) 9418 5333 QLD Laboratory: lab.qld@agrifood.com.au / (07) 3107 9928

